



Spartan Chemical Company, Inc.

1110 Spartan Drive, Maumee OH 43537

Phone: (800) 537-8990

www.spartanchemical.com

Training Information

Market Segment:

Government, Healthcare

Title:

PROCEDURE TO WASH WALLS IN DELIVERY ROOM

Check If Applicable: Rx Procedure None
 Cleancheck Procedure

Tools & Materials Needed:

1. Safety ladder
2. 10-quart bucket
3. Cleaner disinfectant
4. Cleaning cloth
5. Gloves and personal protective equipment

Procedures:

Put on personal protective equipment.

1. Change into scrub suit in doctors' lounge.
2. Fill bucket with detergent solution.
3. Move equipment away from wall to center of Delivery Room.
4. Dip clean cloth in detergent solution and wash ceiling with up and down strokes.
5. Using second bucket and cloth, rinse ceiling. Follow this procedure until all ceiling is washed.
6. Repeat this procedure for the wall, starting at the bottom.
7. Wash overhead lights and tracks at the same time. Using damp cloth, wipe light switches safely and glass partition.
8. Wipe all metal dry with clean cloth to avoid spotting.
9. Mop floor around wall and in all corners. Pick up excess solution.
10. Move equipment back into proper place against wall.
11. Damp mop center of the floor.
12. Repeat process until entire delivery area is completed.